

# From FID to 2D: Processing HSQC Data Using NMRPipe

## Macro Methods

The data recorded during a NMR experiment are stored in the form of a digitized free induction decay (FID), which is in the time domain. In order for us to gain information regarding the chemical environment of individual nuclei that contribute to the observed FID, we need to convert the data into the frequency domain using a Fourier transformation. To do this we will utilize the NMRPipe software package that is maintained by the NIST (<https://www.ibbr.umd.edu/nmrpipe/index.html>). NMRPipe is a powerful program for converting and processing multi-dimensional NMR data. With that in mind, the purpose of this tutorial is to walk you through the various steps involved in processing 2D nitrogen and carbon HSQC datasets recorded for Nb6 and mNb6 alone and in the presence of Spike RBD.

In general, the steps for processing data are: (1) download data to your computer (2) Convert the data to NMRPipe format (3) Fourier transform and phase the first FID and (4) iterate NMRPipe functions through the entire series of 1D experiments.

In order to execute the NMRPipe GUI and commands you will need to either change your default shell to C-shell (csh) or launch a C-shell session in terminal by typing **csh** in a terminal window. In addition, please be sure NMRPipe is installed properly and can be accessed from the terminal—(tutorial for NMRPipe install on MacOS: [http://fraserlab.com/static/pdf/methods/NMRPipe\\_Install\\_MacOS.pdf](http://fraserlab.com/static/pdf/methods/NMRPipe_Install_MacOS.pdf)). To test if things are working, type **nmrPipe** in terminal and it should return the version you have installed.

### Part I: Download data to local computer

1. A Box folder containing all of the data has been shared with you. Download the data into the **nmr** directory you should have created when installing NMRPipe. For ease of locating the data, I would recommend creating a NMR data directory (**mkdir nmrdata**) in your NMR folder.

2. Once the data has been copied over, open terminal (or XQuartz) and descend into the directory:

```
cd ~/nmr/nmrdata/2021_Macro_Methods_NMR
```

**ls** (should now see 5 directories)

```
[Ryans-MacBook-Pro:~/nmr/2021_Macro_Methods_NMR] rwtibble% ls
Nb6                               Processing Scripts               mNb6_RBD_complex
Nb6_RBD_complex                  mNb6
```

We are going to process the HSQC data for Nb6 together and you will then be responsible for processing the remaining datasets in your group.

3. Descend into the directory for Nb6 (**cd Nb6**). Type **ls** again and you should see directories for the carbon (chsqc) and nitrogen (nhsqc) HSQC experiments.

```
[Ryans-MacBook-Pro:~/nmr/2021_Macro_Methods_NMR/Nb6] rwtibble% cd Nb6
[Ryans-MacBook-Pro:~/nmr/2021_Macro_Methods_NMR/Nb6] rwtibble% ls
Nb6_chsqc      Nb6_nhsqc
```

4. We will begin by processing the nitrogen HSQC data (**cd Nb6\_nhsqc**). Type **ls** to look at the contents and ensure the data has copied over correctly:

```
[Ryans-MacBook-Pro:~/nmr/2021_Macro_Methods_NMR/Nb6] rwtibble% cd Nb6_nhsqc/
[Ryans-MacBook-Pro:2021_Macro_Methods_NMR/Nb6/Nb6_nhsqc] rwtibble% ls
EA                format.temp      specpar
acqu              gpnam1          spnam1
acqu2             gpnam2          spnam13
acqu2s           gpnam3          stanprogram2495
acqu              pdata           uxnmr.info
adcInfo_TRX1.xml pulseprogram     uxnmr.par
audita.txt        scon2           vtc_pid_settings
cpdpgr3          ser
format.ased      shimvalues
```


### Part II: Convert from Bruker to NMRPipe format

Next, we need to convert the data from Bruker serial files to NMRPipe format.


5. Type **bruker** in the command line. The NMR Conversion Utility GUI should open.

6. The **Read Parameters** button should be highlighted. Click this to read in the raw data. You should have two columns after clicking this button (x- and y- axis).

# From FID to 2D: Processing HSQC Data Using NMRPipe Macro Methods




7. The parameters marked in the screenshot below, many of which are highlighted in yellow, need to be updated: *Note:* highlighted parameters are obtained from the experimental acquisition parameters and NMRPipe wants you to confirm these are correct.



Click on the arrows by each value highlighted above and select the following options:


- Dimension Count: **2D**
- Center Position PPM (x-axis): **4.700 (O1/BF1) 1H**
- Spectral Width Hz (y-axis): **2127.660**
- Observe Freq MHz (y-axis): **60.803 (SFO3) 15N**
- Center Position PPM (y-axis): **117.00 (O3/BF3) 15N**



# From FID to 2D: Processing HSQC Data Using NMRPipe

## Macro Methods


In addition, an output script we will use for converting the data (titled fid.com) is shown in a second window. Any edits you make in the selection window should also be changed in this script. Make sure the values are correct.




```
#!/bin/csh
bruk2pipe -verb -in ./ser \
-bad 0.0 -ext -aswap -AMX -decim 2080 -dspfvs 21 -grpdly 76 \
-xN 2048 -yN 512 \
-xT 1024 -yT 256 \
-xMODE DDD -yMODE Echo-AntiEcho \
-xSW 9615.385 -ySW 2127.660 \
-xOBS 599.983 -yOBS 60.803 \
-xCAR 4.700 -yCAR 117.000 \
-xLAB HN -yLAB 15N \
-rdim 2 -aq2D Complex \
| nmrPipe -fn MULT -c 7.81250e+00 \
-out ./test.fid -ov

sleep 5
```


8. Click the **Save Script** button (highlighted in green). Click **Continue** on the window that appears to confirm the output file, fid.com, has been written.



9. Next, we will execute the conversion script. Click the **Execute Script** button.




A pop-up window should then appear:



Alternatively, you can run the conversion script by quitting the Conversion Utility and executing the script in the terminal (**`./fid.com`**).

10. NMRPipe should have created a file with the converted data called **test.fid**. Type **ls** and confirm it is present.




```
[Ryans-MacBook-Pro:2021_Macro_Methods_NMR/Nb6/Nb6_nhsqc] rwtibble% ls
EA          format.ased          shimvalues
acqcu      format.temp          spectropar
acqu2     gpnam1              spnam1
acqu2s    gpnam2              spnam13
acqus     gpnam3              stannogram2495
adcInfo_TRX1.xml  pdata                test.fid
audita.txt pulseprogram         uxnmr.info
cpdprg3   scon2                uxnmr.par
fid.com   ser                  vtc_pid_settings
```

# From FID to 2D: Processing HSQC Data Using NMRPipe

## Macro Methods

### Part III: Manually phase the FID using NMRDraw.


11. Open the NMRDraw suite by typing **nmrDraw** in the command line. You should see the following window appear:



*A quick note on navigating NMRDraw:* Use right-click to access the pull-down menu for each of the options in the menu bar. Using left-click automatically performs the first action in each pull down menu. Alternatively, there are letter codes listed to the left of each menu option that can be used to perform a given action.


At this point, NMRDraw automatically read in the test.fid file and displayed the FIDs. You'll notice the scale of the x- and y-axes correspond to the number of real data points we collected for both the hydrogen (direct) and nitrogen (indirect) dimensions.

12. Type **h** to activate the horizontal scroll bar. A pink line should appear along with a single FID.




13. Pull the pink bar to the bottom of the window so that you are visualizing the first recorded FID. The coordinates near the top of the GUI should read **X: 2; Y:1; Z: 1**.

# From FID to 2D: Processing HSQC Data Using NMRPipe Macro Methods




14. Right click the **Proc** button and select **Auto-Process 1D** (or type **p**). This will perform a Fourier transform of the 1-D experiment.




15. Manually phase the spectrum by clicking the phasing **ON** button. Manipulate the zero-order phasing (**P0**) by sliding the left scroll bar back and forth to correct the phasing and make the baseline as flat as possible. The right scroll bar can be used to finely change the phasing.

# From FID to 2D: Processing HSQC Data Using NMRPipe Macro Methods



16. Once you have corrected the phasing, note the number next to **P0** (i.e. 43). Also, make sure **P1** has not been changed from zero.



17. After recording the value for P0, you can exit NMRDraw.

While we just determined the phase correction for a single FID, we need to apply this to all the FIDs we recorded. Luckily, we can run a processing script that will apply the same set of corrections to each FID in order to generate our 2D spectrum.

## Part IV: Implement nhsqc.com file through the data

18. In the Processing\_Scripts directory you downloaded from Box, there is a file called **nhsqc.com**. I recommend keeping 'originals' of this file and other processing scripts in this directory and copy them to your working directory:


```
cp ~/nmr/nmrdata/2021_Macro_Methods_NMR/Processing_Scripts/nhsqc.com .
```

19. Type **ls** to make sure it correctly copied over.

# From FID to 2D: Processing HSQC Data Using NMRPipe Macro Methods

```
[Ryans-MacBook-Pro:2021_Macro_Methods_NMR/Nb6/Nb6_nhsqc] rwtibble% ls
EA                format.temp       specpar
acqu              gpnam1            spnam1
acqu2             gpnam2            spnam13
acqu2s           gpnam3            stanprogram2495
acqu             nhsqc.com         test.fid
adcInfo_TRX1.xml pdata             uxnmr.info
audita.txt       pulseprogram      uxnmr.par
cpdprg3          scon2             vtc_pid_settings
fid.com          ser
format.ased      shimvalues
```

20. Edit the nhsqc.com file using vi or nano (I will use nano here). To do this, type **nano nhsqc.com** and the contents of the file should appear in the terminal window.




```
GNU nano 2.0.6 File: nhsqc.com
#!/bin/csh
#
# Basic 2D Phase-Sensitive Processing:
# Cosine-Bells are used in both dimensions.
# Use of "ZF -auto" doubles size, then rounds to power of 2.
# Use of "FT -auto" chooses correct Transform mode.
# Imaginaries are deleted with "-di" in each dimension.
# Phase corrections should be inserted by hand.

nmrPipe -in test.fid \
nmrPipe -fn SOL
nmrPipe -fn SP -off 0.5 -end 1.00 -pow 1 -c 1.0
nmrPipe -fn ZF -auto
nmrPipe -fn FT -auto
nmrPipe -fn PS -p0 0.00 -p1 0.00 -di -verb
nmrPipe -fn EXT -x1 -left -sw
nmrPipe -fn POLY -auto
nmrPipe -fn TP
nmrPipe -fn SP -off 0.5 -end 1.00 -pow 1 -c 1.0
nmrPipe -fn ZF -auto
nmrPipe -fn FT -auto
nmrPipe -fn PS -p0 -90 -p1 0.00 -di -verb
nmrPipe -fn POLY -auto
-ov -out hsqc.ft2
```

21. Use the arrows to navigate to the PS command and edit the value next to P0 with the one you determined in NMRDraw (i.e. 43)

22. Navigate to the end of the file where it shows the output file (in this case the default is hsqc.ft2) and edit it to describe the spectrum (i.e. Nb6\_nhsqc.ft2). Your edited file should look like this:



```
GNU nano 2.0.6 File: nhsqc.com Modified
#!/bin/csh
#
# Basic 2D Phase-Sensitive Processing:
# Cosine-Bells are used in both dimensions.
# Use of "ZF -auto" doubles size, then rounds to power of 2.
# Use of "FT -auto" chooses correct Transform mode.
# Imaginaries are deleted with "-di" in each dimension.
# Phase corrections should be inserted by hand.

nmrPipe -in test.fid \
nmrPipe -fn SOL
nmrPipe -fn SP -off 0.5 -end 1.00 -pow 1 -c 1.0
nmrPipe -fn ZF -auto
nmrPipe -fn FT -auto
nmrPipe -fn PS -p0 43 -p1 0.00 -di -verb
nmrPipe -fn EXT -x1 -left -sw
nmrPipe -fn POLY -auto
nmrPipe -fn TP
nmrPipe -fn SP -off 0.5 -end 1.00 -pow 1 -c 1.0
nmrPipe -fn ZF -auto
nmrPipe -fn FT -auto
nmrPipe -fn PS -p0 -90 -p1 0.00 -di -verb
nmrPipe -fn POLY -auto
-ov -out Nb6_nhsqc.ft2
```

23. Save this file by pressing **control^X**. Type **Y** to save the modified file and press **ENTER** to overwrite the original nhsqc.com. Alternatively, you can provide a new name and create a new file with this name (becomes more useful as you make more specific changes in the file). If you do choose to create a new file, you will likely have to make it executable by typing **chmod +x nhsqc.com** in the terminal.

24. Run the processing file by typing **./nhsqc.com** in the command line. The file will iterate through all of the experiments and output the ft2 file. If typing this produces an error, you likely need to make the file


## From FID to 2D: Processing HSQC Data Using NMRPipe

### Macro Methods


executable. Do so by typing `chmod +x nhsqc.com` in the terminal and try again. Type `ls` and you should now see the processed 2D file (**Nb6\_nhsqc.ft2**).

```
[Ryans-MacBook-Pro:2021_Macro_Methods_NMR/Nb6/Nb6_nhsqc] rwtibble% ./nhsqc.com
PS      512 of 512   HN
PS      1024 of 1024 15N
[Ryans-MacBook-Pro:2021_Macro_Methods_NMR/Nb6/Nb6_nhsqc] rwtibble% ls
FA      format.ased      shimvalues
Nb6_nhsqc.ft2      format.temp        specpar
acqu    gpnam1             spnam1
acqu2   gpnam2             spnam13
acqu2s  gpnam3             stanprogram2495
acqus   nhsqc.com          test.fid
adcInfo_TRX1.xml  pdata              uxnmr.info
audita.txt         pulseprogram       uxnmr.par
cpdprg3          scon2              vtc_pid_settings
fid.com          ser
```

25. Type `nmrDraw` and you should now see the processed 2D spectrum! If NMRDraw did not automatically open the 2D spectrum, right-click the **File** button and click Select File (or type 'S' from the main screen). A new window should appear. Select `Nb6_nhsqc.ft2` and click **Read/Draw** at the bottom of the window.



26. Use the + and - button in the upper right corner of the window to adjust the display threshold (First: xxxe+xxx) so the peaks more clearly shown. After each click of the + or - button, you'll need to left-click on the Draw button, which resets the contour level (alternatively, type `d`) If you set the contour too low, then you'll start to see many small blue (positive) and red (negative) peaks. Readjust the contour level





## From FID to 2D: Processing HSQC Data Using NMRPipe

### Macro Methods

Look for any indications of a poorly phased spectrum by noticing if peaks appear to be abruptly cut-off and there is a large amount of red, which indicates negative peak height. If this is the case, you can further refine the phasing by typing **h** or **v** to look at individual slices along the x- or y-axes, respectively, and use the scroll bars to improve phasing as we did previously. If you do change the phasing, note this new value, edit, and rerun the processing file. Iterate through this process until you have a well-phased spectrum.